

## Recommended Ordering Information

<b>Sample Collection (choose at least one)</b>	
MicroTest™ M4RT 3mL w/o Beads NP Swab	<a href="#">R12705</a>
<b>Sample Preparation (choose at least one)</b>	
MagMAX™ Viral/Pathogen Nucleic Acid Isolation Kit (Manual extraction, Automated extraction)	<a href="#">A42352</a>
MagMAX™ Viral/Pathogen II (MVP II) Nucleic Acid Isolation Kit (Automated extraction)	<a href="#">A48383R</a>
<b>Automated Nucleic Acid Purification System (optional)</b>	
KingFisher™ Flex Purification System (Saliva automation, VTM automation)	<a href="#">5400630</a>
<b>Manual Nucleic Acid Purification</b>	
PureLink™ Viral RNA/DNA Mini Kit (50 prep)	<a href="#">12280050</a>
MagMAX™ Viral/Pathogen II (MVP II) Nucleic Acid Isolation Kit (2000 prep)	<a href="#">A48383R</a>
<b>Assay (choose SKU based on scale needed)</b>	
TaqMan Monkey Pox Virus Detection Assay (250 rxn)	<a href="#">Vi07922155_s1</a>
<b>Controls</b>	
AcroMetrix Mpox DNA Positive Control RUO	<a href="#">902050</a>
<b>Master Mix</b>	
TaqPath™ BactoPure™ Microbial Detection Master Mix (No ROX™)	<a href="#">A52704</a>
<b>Instruments</b>	
QuantStudio™ 5 Real-Time PCR System, 96-well, 0.2-mL	<a href="#">A28574</a>